

Data sheet

PI / Cell Cycle Analysis Kit

Cat. No: CA112

200 reactions

Introduction

Cell cycle is a ubiquitous, complex sequence of events leading to growth and proliferation of cells. Cell cycle of most eukaryotic cells is divided into four discrete phases: M, G1, S, and G2. S phase is the period during which DNA replication occurs. The cell grows throughout interphase, which includes G1, S, and G2.

Cells at different stages of the cell cycle can also be distinguished by their DNA content. The DNA content of cells duplicates during the S phase (2n to 4n) and DNA content then remains at 4n for cells in G2 and M, decreasing to 2n after cytokinesis.

Experimentally, cellular DNA content can be determined by incubation of cells with a fluorescent dye that binds to DNA, followed by analysis of the fluorescence intensity of individual cells in a flow cytometer. The fluorescence signal is proportional to cellular DNA content.

Propidium iodide is a fluorescent molecule that binds nucleic acid with little or no sequence preference. Because **Propidium iodide** binds RNA as well as DNA, RNase A is included in this kit to digest cellular RNA and thus decrease background RNA staining from the experiment.

Application

- Analysis of cell cycle regulation in response to external stimuli, environmental changes, etc.
- Monitoring of cell cycle progression
- Study effects of drugs which affect cell growth and division.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for Material Safety Data Sheet of the product.

Kit Contents

Components	200 assays
Propidium Iodide (PI) Reagent	2 ml
10x Phosphate Buffered Saline	50 ml
RNase A Solution	2 ml

Storage

Upon receipt, store RNase A solution at -20°C. Store the PI Reagent and 10X PBS at +4°C. Protect PI solution from light. Kit components are stable for at least 12 months from date of receipt when stored as recommended.

Working Solutions

1XPBS: Dilute 10X solutions 10:1 to make a 1x working solution.

Store at room temperature.

Staining Solution: Immediately before use, prepare the Propidium Iodide + RNase A Solution in PBS. For example, for every 20 samples to be stained, mix 10 mL 1XPBS + 200 µL Propidium Iodide + 200 µL RNase A.

Protect from light. Stable for one week at 4°C.

Caution: *Propidium iodide is a suspected carcinogen; contact with eyes, skin, and mucous membranes should be avoided. Always wear proper protective clothing and gloves when handling the solution.*

Canvax Biotech, S.L. C/Astrónoma Cecilia Payne. Edif. Canvax. 14014 Córdoba, Spain.

☎ : +34 957 348 066

☎ : +34 957 346 217

✉ : info@canvaxbiotech.com



Assay Procedure

A. Sample Preparation:

1. Grow cells of interest ($2-5 \times 10^5$ cells/well).
2. Generate a single cell suspension in a manner similar to routine passaging of the cells.
3. Harvest cells and centrifuge at $500 \times g$ for 5 min. Aspirate and discard the supernatant.
4. Wash cells in 1 ml ice cold **1X PBS**.
5. Centrifuge cells at $500 \times g$ for 5 min. Aspirate and discard the supernatant.

B. Nucleic Acid Labelling:


6. Fix the cells by adding 1.2 ml ice cold 70 % ethanol to the cell pellet. Add drop wise to the pellet while vortexing. Put on ice for at least 30 min.

After fixing, cells can be stored at -20°C for several weeks in 70% ethanol

7. Centrifuge cells at $500 \times g$ for 5 min. and carefully remove the supernatant without disrupting the pellet.
8. Wash cells in 1 ml of **1X PBS**.
9. Centrifuge cells at $500 \times g$ for 5 min. and carefully remove the supernatant.
10. Resuspend cells completely with 200 μl of **Staining Solution**, protect from light exposure.
11. Incubate at 37°C for 30 min.
12. Place tubes on ice (still in the dark) and prepare for flow cytometry analysis

C. Data Analysis:

13. Select the main cell population in the FSC vs SSC plot to exclude debris and cell aggregates.
14. Collect propidium iodide fluorescence in the appropriate channel (typically FL2)
15. Setting up markers on a histogram plot to delineate 4N intensity regions facilitates quantifying differences in DNA content between samples within an experiment.



PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for Material Safety Data Sheet of the product.