

Saliva Sample collection Kit

For an easy, convenient and safe collection, transportation and storage of Saliva samples

Ordering info:

Cat No.	Size
SC001	20 units
SC002	100 units
SC003	250 units

Includes for 100 units:

- 100 Funnels
- 100 Saliva Collection Tube
- 100 Saliva Preservation Solution (2.5 mL)



Related products:

- HigherPurity™ Buccal Saliva Genomic DNA Extraction Kit (p.93)
- HigherPurity™ Buccal Swab/Saliva Genomic DNA Extraction Kit (p.93)
- PBS (p.133)
- Custom solutions (p.147)

Description:

The kit offers an easy, convenient and safe collection, transportation and storage of saliva samples. It is based in an optimized Saliva Preservation Solution which allows stabilizes buccal cells and white blood cells found in saliva over 1 year at Room Temperature.

Advantages & Features:

- ✓ **Easy, fast and non-invasive:** sample collection method.
- ✓ **Pure DNA:** ready-to-use in all downstream applications, including sensitive downstream assays.
- ✓ **Convenient:** all-in-one kit, that includes all material needed.
- ✓ **Compatible** with most DNA isolation methods.
- ✓ **Safe:** samples are non-infectious thanks to its optimized Saliva Preservation solution.

Quality control:

- ✓ Rigorous quality control standards to guarantee lot-to-lot consistency.

Applications:

- ✓ PCR.
- ✓ qPCR.
- ✓ Sequencing.
- ✓ SNP analysis.
- ✓ Microarrays.
- ✓ RFLP.
- ✓ Southern Blot Analysis.
- ✓ Methylation studies.