

# Data sheet

## 3 M Sodium Acetate (pH 5.2)

Cat. No BR080: 100 ml

### Description

3M Sodium acetate solution, pH adjusted to 5.2 with glacial acetic acid, Sterile, 100mL.

### Specifications

Chemicals: Biotechnology grade.

Format: Solution.

Concentration: 3M .

Volume: 100 ml.

pH: 5.2 ± 0.05 at 25 °C.

### Shipping and storage

Sodium Acetate is shipped at room temperature. Store the product at Room Temperature.

### Directions for use

Ready to use solution.

### Applications

Ideal for a wide range of molecular biology applications.

Canvax Biotech, S.L. C/Astrónoma Cecilia Payne. Edif. Canvax. 14014 Córdoba, Spain.

☎ : +34 957 348 066

☎ : +34 957 346 217

✉ : [info@canvaxbiotech.com](mailto:info@canvaxbiotech.com)



Research Use Only (RUO)